

Kinser, Robin D.

From: Kinser, Robin D.
Sent: Monday, October 07, 2002 6:50 PM
To: Miert, Erik Van; Vanscheeuwijck, Patrick
Cc: Kinser, Robin D.; Feng, Shixia; Roethig, Hans
Subject: 3-HPMA collection

Dear Erik and Patrick--

As promised, I have spoken with the analyst who developed the 3-HPMA analysis at the contract research laboratory we employ. He agrees that the brief exposure of the urine to the untreated glass during the voiding time of the rat should not cause a problem. He does encourage you to chill the collected urine as quickly as possible--hopefully within 3 hours. They have seen a 15% decrease in samples after 6 hours of room temperature exposure.

I will be on holiday until October 21; if you have further questions you may contact my colleague Shixia Feng at 804-274-3758 or Shixia.Feng@pmusa.com.

Best regards,
Robin

Robin Kinser
804-274-5905
Robin.D.Kinser@pmusa.com